BORT ManuStabil kurz

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REF 112 730

D112730|2020-10|004 ML

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DE BORT ManuStabil kurz

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EN BORT ManuStabil short

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose
This medical device is a wrist brace for immobilisation in at least two directions of movement.

Indications
Pre-operative, post-operative, post-traumatic, e.g. in case of distortion, degenerative complaints, e.g. carpal tunnel syndrome (CTS), Tenosynovitis, peripheral nerve paralysis.

Contraindications
Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

Application risks/Important notes
This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
– remove the medical device prior to radiological examinations
– do not close it too tightly during longer periods of rest, e.g. when sleeping
– in case of numbness, loosen the medical device and remove it if necessary
– in case of persistent complaints, consult the physician or specialist retailer
– use the medical device in accordance with therapeutic needs
– only use other products simultaneously after consultation with your physician
– do not make any changes to the product
– do not wear it on open wounds
– do not use in case of intolerance of one of the materials used
– no re-use – this medical aid is intended for treating one patient
– whilst wearing the brace, please neither use any creams nor ointments in the vicinity of the attached aid as they can destroy the material

Fitting instructions
Before putting on, bend the aluminium splints to the desired shape. Position the dorsal splint in one of the two available splint bags, dependent on the indication position. Open the Velcro fasteners and the strap and place the brace on the wrist affected. Close the fastener between the thumb and index finger first, then the following Velcro fasteners. Wind the strap circumferentially around the wrist and attach it to the velour using the Velcro.

Removal
To remove the brace, please open the strap and the Velcro fasteners. Grasp

the splint in the palm area and pull the brace off over the fingers. Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

Material composition
black: polyamide (PA), polyester (PES), cotton (CO), elastane (EL)
grey: polyamide (PA), polyester (PES), elastane (EL)
skin-tone: polyester (PES), polyamide (PA), viscose (CV), elastane (EL), polyurethane (PUR)
The sewn in textile label provides the precise material composition.

Cleaning information
 Delicate wash Do not clean chemically Do not bleach
 Do not dry in a tumble dryer Do not iron
Remove the aluminium splints before washing. Close the Velcro fasteners in order to avoid damaging other items of clothing. Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

Guarantee
The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product
The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification
If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA). You can find our contact information in these instructions for use. You can find the contact information for the appointed authority for your country under the following link: www.bort.com/md-eu-contact.

Disposal
Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity
We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 11.2019

Medical device | Single patient – multiple use

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